


General Protocol for Diluted* Intervention™ Concentrate:

1

Clean (remove organic matter)



3

Apply disinfectant (foaming recommended)




2

Rinse if using a detergent



4

Wait 5 minutes to be confident disinfection has occurred. No need to rinse.




* In sub-freezing temperatures, replace up to 40% of water volume with antifreeze (propylene glycol or PG) without affecting AHP's efficacy

General Protocol for Intervention™ Ready to Use (RTU) Liquid:

1

Clean (remove organic matter by using a brush, shaking out floor mats, or washing with detergent)



2

Apply using a trigger foamer



3

Wait 60 seconds to be confident disinfection has occurred

